


ION TEST Analysis Procedure



△Read carefully Manual in advance
△Keep ION TEST Kits in a cool and dark place

Where and Which Parameter is IONTEST applied?


for Environmental Research

COD, NO₂, NO₃, NH₄, PO₄, TN-i, etc.




for Drinking Water Inspection

O₃, HOCl, COD, NO₂, NO₃, NH₄, TN-i, Fe, etc.



for Plating Analysis

Copper, Nickel, Sulfate plating solution industry.
Cu, Ni, Cr⁶⁺, CN⁻, PO₄, Fe, Zn, H₂O₂, etc.



for Waste Water Analysis

Food Factory, Chemical Factory, Drainage, etc.
COD, TN-i, PO₄, NO₂, NO₃, NH₄, HOCl, O₃, Cr⁶, Cr³, CN⁻, Zn, etc.



for Fishery Management

COD, NH₄, NO₂, NO₃, PO₄, TN-i, etc.



for Educational Aid

COD, NO₂, NO₃, NH₄, TN-i, Cu, Ni, Cr⁶, HOCl, O₃, CN, H₂O₂, Fe, Zn, etc.



KRK

- ◆ Powder Reagent is sealed up for long /safe validity
- ◆ Convenient for on-site Testing
- ◆ Refined Colorimetry, Ion-Test
Fine transparent Test Cell & 7 colors gradation


Colorimetric ION TEST

Water Simple Analysis Kit



KRK **Kasahara**
CHEMICAL INSTRUMENTS CORP.
2-133-8 Sakurada Kuki-City Saitama Japan 〒340-0203
TEL +81-480-38-9151 FAX +81-480-38-9157
www.krkjpn.co.jp
E-mail.krk@krkjpn.co.jp

Kasahara



IONTEST

Water Simple Analysis kit

Anyone can Easily Analyze Water Quality, Anytime and Anywhere.

ION TEST

| Model | Testing item | Range (mg/ℓ = ppm) | Determination Time | Quantity | Applications |
|-------------------|-----------------------------------|--|--------------------|----------|--|
| WIT-Cu | Copper | 0.2/0.5/1.0/2.0/3.0/5.0/10 Reduction and Bicinchoninic Acid method | 2min. | 50 | Waste water management |
| -Cu(B) | Copper | 0.2/0.5/1.0/2.0/3.0/5.0/10 Bathocuproine Disulfonic Acid Disodium Salt method ★ If Test water contains such Chelating agent as EDTA, etc., please apply WIT-Cu(B) for testing. | 1min. | 50 | Waste water management |
| -Ni | Nickel | 0.2/0.5/1.0/2.0/3.0/5.0/10 Nioxime method | 2min. | 50 | Waste water management |
| -NH ₄ | Ammonium ion Ammonium-Nitrogen | 0.3/0.7/1.3/2.6/6.5/13/26 NH ₄ ⁺ 0.2/0.5/1.0/2.0/5.0/10/20 NH ₄ ⁺ -N Indophenol Blue method | 5min. | 50 | Waste water management River water check |
| -COD-M | COD-M(Middle range) | 0/5/10/13/20/50/100 Oxidation with Potassium Permanganate in Alkalinity method | 4~6min. | 50 | Waste water management River water check |
| -COD-H | COD-H(High range) | 0/20/40/80/120/180/250 Oxidation with Potassium Permanganate in Alkalinity method | 4~6min. | 50 | Waste water management Environmental research |
| -Cr ⁶⁺ | Chromium(Hexavalent) | 0.05/0.1/0.2/0.5/0.8/1.0/2.0 Cr ⁶⁺ Diphenylcarbazide method | 2min. | 50 | Waste water management Process management |
| -Cr ^T | Total Chromium | 0.5/1/2/5/8/10/20 Cr ^T Oxidation and Diphenylcarbazide method | 30sec. | 50 | Waste water management RoHS |
| -HOCl | Residual Chlorine | 10/20/30/50/80/100/150 KI method | 10sec. | 50 | Waste water management Residue check |

| Model | Testing item | Range (mg/ℓ = ppm) | Determination Time | Quantity | Applications |
|------------------------------------|--|---|--------------------|----------|--|
| -O ₃ | Ozone | 0.1/0.2/0.3/0.5/0.8/1.0/2.0 DPD method | 10sec. | 50 | Waste water management |
| -CN | Free Cyanide | 0.02/0.05/0.1/0.2/0.5/1.0/2.0 4-Pyridinecarboxylic Acid method | 10min. | 50 | Waste water management Poison detection |
| -H ₂ O ₂ -H | Hydrogen Peroxide (High range) | 10/20/30/50/80/100/150 KI method | 10sec. | 50 | Waste water management |
| -NO ₂ | Nitrite Nitrite-Nitrogen | 0.05/0.1/0.2/0.4/0.6/0.8/1.0 NO ₂ ⁻ 0.015/0.03/0.06/0.12/0.18/0.24/0.30 NO ₂ ⁻ -N GR method | 3min. | 50 | Waste water management River water check |
| -NO ₃ | Nitrate Nitrate-Nitrogen | 0.5/1.0/2.0/4.0/6.0/10/20 NO ₃ ⁻ 0.1/0.2/0.5/1.0/1.4/2.3/4.6 NO ₃ ⁻ -N Reduction GR method | 3min. | 50 | Waste water management River water check |
| -PO ₄ | Phosphate Phosphate-Phosphorus | 0.2/0.5/1.0/1.5/2.0/3.0/5.0 PO ₄ ³⁻ 0.1/0.2/0.3/0.5/0.7/1.0/1.7 PO ₄ ³⁻ -P Molybdenum blue method (Measurement with 2 reagents) | 3min. | 50 | Waste water management River water check |
| NEW -PO ₄ (B) | Phosphate Phosphate-Phosphorus | 0.6/1.5/3.0/4.5/6.0/9.0/15 PO ₄ ³⁻ 0.2/0.5/1.0/1.5/2.0/3.0/5.0 PO ₄ ³⁻ -P Molybdenum blue method (Measurement with 1 reagent) | 3min. | 50 | Waste water management River water check |
| -PO ₄ -H | Phosphate Phosphate-Phosphorus (High range) | 2/5/10/15/20/30/50 PO ₄ ³⁻ 0.7/1.7/3.3/5.0/6.6/10/17 PO ₄ ³⁻ -P Molybdenum blue method | 3min. | 50 | Process management Environmental research |
| -TN-i | Total Nitrogen (Inorganic) | 0/5/10/20/40/60/100 Reduction & Indophenol blue method | 20min. | 50 | Waste water management River water check |
| -Fe | Iron | 0.2/0.5/1.0/1.5/2.0/3.0/5.0 Reduction & O-phenanthroline method | 5min. | 50 | Waste water management Drinking water check |
| -Zn | Zinc | 0/0.2/0.3/0.5/1.0/2.0/5.0 PAN method | 2min. | 50 | Waste water management |

ION TEST KIT SPECIFICATION

Standard Component

| | | |
|---|----------------|---------|
| 1 | Standard Color | 1 |
| 2 | Reagent | 50 Pcs. |
| 3 | Test Cell | 1Pc. |
| 4 | Box | 1Pc. |
| 5 | Pipette | 1Pc. |



STD Component part is optionally available.

| | | |
|---------------------------------|------------------------------|--------------------------|
| | | |
| Standard Color (selective) | Reagent (50 pcs) (selective) | Test Cell with Lid (5ml) |
| WIT-□□□□-STD (EXP. WIT-Cu-STD.) | WIT-□□□□ (EXP. WIT-Cu) | CELL-5C |

Features

- Test Cell/Color Card of repeatable use
- Powdered reagent is sealed air-tight & economically priced.
- Colorimetry in 7 colors
- Anyone can easily analyze anytime, anywhere

Test Application

- Semiconductor, Plating industry sewage
- Tap water, river water, groundwater
- Experimental chemistry at laboratory, school.
- Septic tank drainage, Factory drainage